14307 Vitro Fil LC - Instrucao de Uso

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VITRO FIL LC
Glass Ionomer Cement
English

PROFESSIONAL USE ONLY

VITRO FIL LC is a triple, resin-modified cement of glass ionomer of easy handling to be used in direct restorations. It is light-activated by visible light which facilitates the operation. Light-activation provides a longer working time and a fast jellification. The light-activation allows the rapid stabilization of the material (20 seconds), minimizing the sensitivity to water during the initial stages of jellification.

VITRO FIL LC offers a strong bonding of the ionomer to the tooth. It releases fluoride in the dental structure and the fluoride absorption minimizes the incidence of secondary caries. Its excellent radiopacity facilitates the radiographic control. The cement of glass ionomer offers high chemical adhesion, fluoride release, biocompatibility and good esthetic.

Composition:

VITRO FIL LC powder contains fluorine strontium aluminum silicate, excipients, activators and iron oxide.

VITRO FIL LC liquid contains 2-hydroxyethyl methacrylate, polyacrylic and tartaric acid solutions, benzoyl peroxide and camphorquinone.

VITRO FIL PRIMER contains modified methacrylate polyacids, stabilizer, catalyst and ethyl alcohol.

NATURAL GLAZE contains bisphenol glycidyl methacrylate, triethylene glycol dimethacrylate, 2,6-terc-butylphenol, ethyl urethane, B200P, benzyl dimethyl ketal, camphorquinone and quantacure EHA.

Indications:

It is indicated for Class III and V restorations and particularly in restorations of cervical erosions and of root surfaces caries. For restorations in deciduous teeth, as a liner and as filling cores (mainly in vital teeth). It can be also used in small Class I and II restorations, where there are not occlusal contact points.

Contraindications:

It must not be used over exposed pulp or on dentine which pulp thickness is less than 0.5 mm. Avoid the use of VITRO FIL LC in patients with previous history of allergy to metacrylates.

Presentations:

Kit with 5g bottle of VITRO FIL LC powder + 2.5mL bottle of VITRO FIL LC liquid + 2.5mL bottle of VITRO FIL PRIMER + 5mL bottle of NATURAL GLAZE + 1 measuring spoon + 1 mixing pad. Available in the following shades: A1, A2, A3, A3.5, B2, B3 and C3.

Mini kit with 5g bottle of VITRO FIL LC powder + 2.5mL bottle of VITRO FIL LC liquid + 1 measuring spoon + 1 mixing pad. Available in the following shades: A1, A2, A3, A3.5, B2, B3 and C3.

Box with a 5g bottle of VITRO FIL LC powder + 1 measuring spoon. Available in the following shades: A1, A2, A3, A3.5, B2, B3 and C3.

Box with a 2.5mL bottle of VITRO FIL LC liquid.

Storage and Conservation:

VITRO FIL LC powder and liquid, VITRO FIL PRIMER and NATURAL GLAZE must be stored between 15°C and 30°C. Do not expose to excessive heat. Keep the containers tightly closed after use. Under these conditions, the product is valid for up to 2 years as of the manufacturing date. Do not use this product if shelf-life has already expired.

Instructions for Use:

The use of gloves is recommended during the procedure. Attention: The patient must be informed about the risk of sensitivity in the place of the procedure. This sensitization can be minimized through the use of proper isolation (rubber dam) and pulp protection.

Cavity Preparation

Perform the usual cavity preparation.

Cavity Cleaning and Drying

The site that will receive the preparation should be cleaned with pumice and water. Rinse it thoroughly and dry. Avoid dehydration.

NOTE: The use of VITRO CONDICIONADOR (CONDITIONER) for 10 seconds is recommended to remove the dentine sludge. Rinse and dry, but make sure to avoid dehydration. Apply calcium hydroxide in areas of pulp exposure or dentine with thickness less than 0.5mm. To achieve extra adhesion, as in the case of a Class V restoration, you should apply VITRO FIL PRIMER and light-activate for 20 seconds. On this type of union, the adhesion occurs by means of a micromechanics bonding with dentine (dentinary adhesion of 3rd generation) and not by means of Van der Waals strength, proper of ionomeric adhesion.

Dosage

The proportion used should be of 1 measuring spoon of the powder to 2 drops of the liquid. Before take the powder, the bottle should be gently shaken. Fill in the measuring spoon with powder and level by taking the excess. The liquid should be dispensed holding the bottle in a vertical position, with the neck turned down. Any crystals present in the neck should be removed to assure that there are not air bubbles in the drop dispensed. Keep the powder and liquid bottles tightly closed after use.

Mixture

Add the necessary amounts of powder and liquid to the mixing pad. Using the plastic or metal spatula, divide the powder into two equal portions. Join the liquid to the first portion of powder and mix for 10 seconds. After that, add the other portion of powder and mix for 10-15 seconds. The total mixture time should not exceed 20-25 seconds. The material is indicated to be used in a temperature between 21-30°C. The handling time is of 3 minutes after the finish of the mixture at 23°C. Higher temperatures will reduce the working time.

Application and Initial Carving

a) Use different mixtures for each element to be restored. With a spatula #1 or a Holleback take

a small amount of the mixed material and apply the cement on the preparation site avoiding air bubbles.

Another option to avoid air bubbles is use material dispenser with the corresponding tip.
b) After applying the cement, contours and sculpture should be performed. A polyester band or a matrix over the freshly applied cement can be used. It can be removed after material jellification.

Light-activation

Light cure for 20 seconds using light-curing devices with a minimum intensity of 820mW/cm². Devices with lower light intensity may require a longer curing time.
NOTE: If the application site exceeds 2mm depth, the lining technique should be used.

Finishing

a) After cement hardening, finishing can be accomplished with a water jet, if necessary. Use an extra-thin diamond tip or silicone tip.
b) After finishing apply NATURAL GLAZE, light curing for 20 seconds to seal restoration surface.

Properties of Use

Times given below are valid for the specific temperature and humidity conditions described. Higher temperatures will reduce and lower temperatures will increase working time.

Time of work (Low ambient light, 21-23°C, 50-52% HR)	> 2 min
Time of work (In the dark, 21-23°C, 50-52% HR)	4 min - 7 min
Time of cure (36-38°C, > 90% HR)	2 min 30 sec - 5 min

The product tends to adhere to metal instruments; for this reason, it is recommended that these instruments be cleaned with cold water before residues drying. Product wastes must not be reused. They must be disposed as the local law in force.

Precautions:

Avoid the use of VITRO FIL LC in patients with previous history of allergy to glass ionomer cements. In case of allergy, you must discontinue the use immediately and advise the patient to see a doctor. The liquid and the prepared cement can be corrosive. Do not allow contact with soft tissues or the skin. In case of contact, you should immediately remove the material with hydrophilic cotton soaked with water and rinse thoroughly with water. Avoid contact with the eyes; in case of accidental contact, wash the eyes immediately with water and see a doctor. Do not mix the powder or liquid which are components of this product with any other glass ionomer-based product. This product must be used only by surgeon dentists and for the intended use described in the Instructions for Use. Avoid contact with eugenol-based products, since they can delay the action of glass ionomers.

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DFL Indústria e Comércio S.A.