150mm

CE
2460

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VITRO CEM
Glass Ionomer Luting Cement
English

PROFESSIONAL USE ONLY

Composition:

The powder contains: Fluorine Strontium Aluminium Silicate, Dehydrated Polyacrylic Acid and Iron Oxide.
The liquid contains: Polyacrylic Acid, Tartaric Acid and Distilled Water.
VITRO CEM is a chemically activated acid/base glass ionomer cement. It is radiopaque due to the strontium action and is easily manipulated. It presents the features that are common to other glass ionomers, i. e., biocompatibility, high chemical adhesion, low solubility, adhesion to wet surfaces and continuous fluoride release.

Indications:

VITRO CEM is recommended for cementation of inlays/onlays, crowns, bridges, prefabricated pins, melted nuclei, orthodontic bands and for the use as liner or base under dental restorations. It presents high radiopacity and fast hardening.

Presentations:

- Box containing: 1 bottle of **VITRO CEM** powder (5g) + 1 bottle of **VITRO CEM** liquid (4mL) + 1 measuring spoon.
- Box containing: 1 bottle of **VITRO CEM** powder (10g) + 1 bottle of **VITRO CEM** liquid (8mL) + 1 measuring spoon + 1 mixing pad.
- Box containing: 1 bottle of **VITRO CEM** powder (15g) + 1 bottle of **VITRO CEM** liquid (10mL) + 1 measuring spoon + 1 mixing pad.
- Box containing: 1 bottle of **VITRO CEM** powder (10g) + 1 measuring spoon.
- Box containing: 1 bottle of **VITRO CEM** powder (15g) + 1 measuring spoon.
- Box containing: 1 bottle of **VITRO CEM** powder (33g) + 1 measuring spoon.
- Box containing: 1 bottle of **VITRO CEM** liquid (8 mL).
- Box containing: 1 bottle of **VITRO CEM** liquid (10 mL).

Storage and Conservation:

VITRO CEM must be stored between 15°C and 30°C. Keep the bottles of powder and liquid tightly closed after use. Do not expose to excessive heat. Under these conditions, the product is valid for up to 2 years as of the manufacturing date. Do not use this product if shelf life had expired.

Instructions for Use:

The use of gloves is recommended during the procedure.
Attention: The patient must be informed on the risk to occur sensitivity in the place of the procedure. This sensitization can be minimized through the use of proper isolation (rubber sheet) and pulp protection.

To obtain maximum adhesion, it is recommended to clean and dry the enamel, dentin and metal surfaces. However, excessive drying of dental structures must be avoided. A wet pellicle, as frequently used in vital teeth, can be incorporated to the cement.

Dosage

Gently shake the bottle, to release the powder. Fill the measuring spoon with powder and level it in the stopper. The proportion used should be of 1 measuring spoon of the powder to 2 drops of the liquid. The liquid should be dispensed while holding the bottle in a vertical position. Any crystals present in the neck should be removed, to assure that there are not air bubbles in the dispensed drop. Keep the powder and liquid bottles tightly closed after use, avoiding exposure to light.

Mixture

With a metal or plastic spatula, prepare the mixture over an appropriate mixing pad or glass slab during 20 seconds. The ideal procedure is to add the powder into the liquid in 2 portions, adding the second portion just after complete homogenization of the first portion. Mix the powder with the liquid until a homogenous consistence is obtained. For the collage of a crown, the amount of cement should be taken in the spatula in a single time. **VITRO CEM** must be manipulated at a temperature between 20-30°C. Times given below are valid for the specific temperature and humidity conditions described. Higher temperatures will **reduce** and lower temperatures will **increase** working time.

Mixing Time (21-23°C, 50-52% UR)	20 sec
Working Time (21-23°C, 50-52% UR)	> 2 min 30 sec
Net Setting Time (36-38°C, = 90% UR)	3 to 4 min

Attention: A higher powder/liquid ratio also reduces the working time. If a longer working time is desired, cool the glass slab. An excessively thin mixture (due to little quantity of powder) can cause a higher solubility and therefore, lead to hyperesthesia symptoms.

Application

During application and curing, you should avoid contact with water or saliva. It is recommend to apply calcium hydroxide immediately after preparation, in cases where the residual dentin layer is too thin and/or in the presence of exposed pulp. Avoid excessive materials in crowns. Brush a thin layer of cement on the internal surface of the crown. Place the prosthesis in the definite position, without applying any pressure. Hydrostatic pressure in excessively tight crowns can cause pulp sensitivity. If working time is exceeded, cement adhesion to enamel and dentin can be lost. The product tends to adhere to metal instruments; for this reason, it is recommended that these instruments be cleaned with cold water before residues drying.
Product wastes must not be reused. They must be disposed as the local law in force.

Precautions

The liquid and the prepared cement can be corrosive. Do not allow contact with soft tissues or the skin. In case of contact, you should immediately remove the material with hydrophilic cotton soaked with water and rinse thoroughly with water. Avoid contact with the eyes; in case of accidental contact, wash the eyes immediately with water and see a doctor. Do not mix the powder or liquid which are components of this product with any other glass ionomer-based product. This product must be used only by dentists and for the intended use described in the Instructions for Use. Avoid contact with eugenol-based products, since they can delay the action of glass ionomers.

Contraindications:

Avoid the use in patients with previous history of allergy to glass ionomer cements. In case of allergic reaction, discontinue the use immediately and advise the patient to see a doctor.

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DFL Indústria e Comércio S.A.